

# Explore

## Ultra transportable power wheelchair



*Available Spring 2012*

The Explore is a compact and lightweight power wheelchair. It can be disassembled in one minute (without tools) into three modules, to fit into the trunk or backseat of any vehicle. No need for an expensive van! It can also be packed in an optional suitcase for travel.

Its compact design means that the Explore is extremely manoeuvrable in tight spaces and all of this is offered with standard wheelchair functionality.

### Specifications:

- **Weight:** 24 kg. The front frame with battery weighs 6 kg. The rear frame, with motor weighs 11 kg. The seat and controller weigh 6 kg.
- **Battery:** The 24 volt 10 AH Lithium Polymer battery weighs 1.8 kg. Its range is 19.2 km (on flat terrain).
- **Seat:** The standard seat width and depth are 16 inches (406x415 mm). 18"x16" (456x415 mm), 14"x15" (352x383 mm), and 12"x15" (302x383 mm) seat sizes are also available. The seat backrest angle can be set between 90 to 105 degrees. The seat angle can be adjusted from 0-5 degrees.
- **Wheels:** The 9 inch (2.8/2.50-4) tires and suspension give the Explore outstanding drivability
- **Personalization:** The foot rest height, position and angle, as well as the arm rest height, can be adjusted according to the height of the user. The wheelbase length can be set to between 430 mm and 530 mm for extra stability.
- **Speed:** 8.4 kph
- **Weight capacity:** 115 kg.

### Contact us :

Henry Dunantstraat 45 | 7548 AB Enschede | Netherlands  
Tel. : + 31 (0)53-4281589  
Mob. : + 31 (0)640-596959  
Fax : + 31 (0)53-4500453  
Internet: [www.rollick.biz](http://www.rollick.biz) | [info@rollick.biz](mailto:info@rollick.biz)

